

SD-B
Injector nozzle
for butyrometers.



- For mobile injector units
- Including thread, 4 x 140 mm*

* Including thread, diameter 4 x 140 mm, plus welded flattened nozzle, 1.5 x 100 mm.

EAN: 4002511116955 / Material number: 03583540 / Old Material Number: 69708001D

Product classification	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category	
Injector nozzle	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	240
External dimensions, gross height in mm	12
External dimensions, gross width in mm	100
External dimensions, gross depth in mm	150
Diameter in mm	4
Net weight in kg	0.01
Gross weight in kg	0.02