


O 183

Upper basket

for the optimum loading of instruments for MIS, arthroscopy and urology.



- Left side: 10 injector nozzles or free accommodation of inserts
- Right side: 14 injector nozzles/ funnels for hollow instruments
- Vertical clearance 285 +/- 30mm; hollow instruments 370 +/- 30mm
- Powder dispensing not possible

EAN: 4002511114920 / Material number: 03807470 / Old Material Number: 69118301D

| Product category | |
|---|-----------------|
| Washer-disinfectors, medical | • |
| Machine affiliation | |
| PG 8535 | • |
| PG 8536 | • |
| Product category | |
| Injector basket | • |
| Upper basket | • |
| Product properties | |
| Material | Stainless steel |
| Colour | Stainless Steel |
| Area of application | |
| Reprocessing of MIS instruments | • |
| Reprocessing of lumened instruments | • |
| Capacity | |
| Lumened instruments [number] | 14 |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0 |
| Frequency in Hz | 0 |
| Total rated load in kW | 0.00 |
| Fuse rating in A | 0 |
| Dimensions and weight | |
| External dimensions, net height in mm | 270 |
| External dimensions, net width in mm | 475 |
| External dimensions, net depth in mm | 530 |
| External dimensions, gross height in mm | 410 |
| External dimensions, gross width in mm | 540 |
| External dimensions, gross depth in mm | 555 |
| Net weight in kg | 4.00 |
| Gross weight in kg | 4.80 |