


E 106

Insert

for the optimum loading of wide neck flasks, measuring cylinders, etc.



- Equipped with 10 spring hooks x 175 mm height and 16 x 105 mm height
- The distance between the spring hooks is approx. 60 mm

EAN: 4002511115675 / Material number: 03808310 / Old Material Number: 69510601D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 8505	•
Product category	
Insert/module for beakers	•
Insert/module for wide-necked glassware	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Reprocessing of measuring cylinders	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	180
External dimensions, net width in mm	195
External dimensions, net depth in mm	186
External dimensions, gross height in mm	205
External dimensions, gross width in mm	245
External dimensions, gross depth in mm	490
Net weight in kg	0.93
Gross weight in kg	1.42