

E 106/1

Insert

For the optimum loading of wide neck flasks, measuring cylinders, etc.



- Equipped with 26 small spring hooks x 105 mm height
- The distance between the spring hooks is approx. 60 mm

EAN: 4002511115682 / Material number: 03808320 / Old Material Number: 69510602D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 8505	•
<b>Product category</b>	
Insert/module for beakers	•
Insert/module for wide-necked glassware	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
Reprocessing of measuring cylinders	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	195
External dimensions, net width in mm	116
External dimensions, net depth in mm	410
External dimensions, gross height in mm	250
External dimensions, gross width in mm	190
External dimensions, gross depth in mm	490
Net weight in kg	0.77
Gross weight in kg	0.94