


E 109

Insert

for the optimum loading of beakers up to 250 ml.



- Suitable for 21 glass beakers, max. 250 ml
- equipped with 21 x 3 stay bars
- For use in upper and lower basket

EAN: 4002511115729 / Material number: 03808360 / Old Material Number: 69510901D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 8505	•
<b>Product category</b>	
Insert/module for beakers	•
<b>Product properties</b>	
Material	Stainless steel, Wire, rilsan-coated
Colour	Stainless Steel, White
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Beakers [quantity]	21
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	155
External dimensions, net width in mm	230
External dimensions, net depth in mm	460
External dimensions, gross height in mm	185
External dimensions, gross width in mm	230
External dimensions, gross depth in mm	490
Net weight in kg	1.21
Gross weight in kg	1.54