

E 129

Insert

For the optimum loading of infusion flasks. Total area 84 x 84 mm.



- 20 compartments, dimension: 84 x 84 mm (bottle), 46 x 46 mm (neck)

EAN: 4002511115866 / Material number: 03808550 / Old Material Number: 69512901D

<b>Product classification</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
PG 8562	•
<b>Product category</b>	
Insert for various flasks	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	20
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	113
External dimensions, net width in mm	445
External dimensions, net depth in mm	445
External dimensions, gross height in mm	480
External dimensions, gross width in mm	480
External dimensions, gross depth in mm	130
Net weight in kg	1.41
Gross weight in kg	1.78