

E 135/1

Insert

for optimum loading of 19 baby bottles up to 110 ml.



- Suitable for 19 baby bottles x 110 ml
- Suitable for bottles sized 51 x 51 mm, neck: 45 x 45 mm
- For upper basket A 102 and lower basket A 151

EAN: 4002511115910 / Material number: 03808640 / Old Material Number: 69513502D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7826	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Insert for various flasks	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of baby bottles	•
<b>Capacity</b>	
Baby bottles [number]	19
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	135
External dimensions, net width in mm	192
External dimensions, net depth in mm	447
External dimensions, gross height in mm	195
External dimensions, gross width in mm	205
External dimensions, gross depth in mm	465
Net weight in kg	2.02
Gross weight in kg	2.10