


E 149

Insert

for optimum loading of 80 test tubes, 16 x 105 mm.



- Suitable for 80 test tubes up to 16 x 105 mm
- Additionally equipped with a lid
- 80 compartments at 18 mm x 18 mm
- Mesh width on the base 8 x 8 mm

EAN: 4002511116016 / Material number: 03808800 / Old Material Number: 69514901D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 8505	•
<b>Product category</b>	
Insert/module for test tubes	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	80
Compartment width in mm	18
Compartment length in mm	18
<b>Area of application</b>	
Reprocessing of test tubes	•
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Test tubes [number]	80
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	152
External dimensions, net width in mm	200
External dimensions, net depth in mm	200
External dimensions, gross height in mm	170
External dimensions, gross width in mm	250

E 149  
Insert  
for optimum loading of 80 test tubes, 16 x 105 mm.

EAN: 4002511116016 / Material number: 03808800 / Old Material Number: 69514901D

External dimensions, gross depth in mm	250
Net weight in kg	1.03
Gross weight in kg	1.59

