

E 353

Clamps

Height-adjustable for optimum positioning of laboratory glassware on the nozzle.

EAN: 4002511117099 / Material number: 03809530 / Old Material Number: 69735301D



- Suitable for injector nozzle E 351
- Height adjustable locking
- Suitable for bottles, Erlenmeyer flasks and round bottomed flasks for example

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category	
Other accessories	•
Product properties	
Material	Wire, rilsan-coated, Stainless steel
Colour	Stainless Steel, White
Area of application	
Reprocessing of laboratory glassware	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	100
External dimensions, net width in mm	45
External dimensions, net depth in mm	160
External dimensions, gross height in mm	120
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	100
Net weight in kg	0.03
Gross weight in kg	0.04