

E 354

Clamps

Height-adjustable for optimum positioning of laboratory glassware on the nozzle.

EAN: 4002511117105 / Material number: 03809540 / Old Material Number: 69735401D



- Suitable for injector nozzle E 352
- Height adjustable locking
- Suitable for bottles, Erlenmeyer flasks and round bottomed flasks for example

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
<b>Product category</b>	
Other accessories	•
<b>Product properties</b>	
Material	Wire, rilsan-coated, Stainless steel
Colour	Stainless Steel, White
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	150
External dimensions, net width in mm	70
External dimensions, net depth in mm	70
External dimensions, gross height in mm	160
External dimensions, gross width in mm	100
External dimensions, gross depth in mm	150
Net weight in kg	0.05
Gross weight in kg	0.05