

E 374

Insert

for optimum loading of ENT instruments such as inhalation connections.



- 24 comp. approx. 45 x 45 mm, 27 compartments approx. 12 x 12 mm
- Mesh width on base: 1.7 mm
- For use in upper and lower basket

EAN: 4002511116412 / Material number: 03809800 / Old Material Number: 69537401D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Insert for ENT instruments	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	51
<b>Area of application</b>	
Reprocessing of ENT instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	63
External dimensions, net width in mm	173
External dimensions, net depth in mm	445
External dimensions, gross height in mm	75
External dimensions, gross width in mm	190
External dimensions, gross depth in mm	460
Net weight in kg	1.20
Gross weight in kg	1.20