

E 408  
Mobile unit  
for the optimum loading of up to 96 pipettes.



- Suitable for 96 pipettes up to 450 mm
- Surface area 16 x 16 mm
- Complete with frame height 150 mm
- Equipped with connection for hot air drying

EAN: 4002511116689 / Material number: 03810130 / Old Material Number: 69540801D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
<b>Product category</b>	
Mobile unit with connection for hot-air drying	•
Mobile unit for pipettes	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Area of application</b>	
Reprocessing of pipettes/phials	•
<b>Capacity</b>	
Pipettes [number]	96
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	502
External dimensions, net width in mm	535
External dimensions, net depth in mm	517
External dimensions, gross height in mm	540
External dimensions, gross width in mm	530
External dimensions, gross depth in mm	500
Net weight in kg	5.50
Gross weight in kg	7.10