


E 413

Insert

for the optimum loading of extraction hoses.



- For 6 System Sirona suction tubes
- For use in upper basket O177/1

EAN: 4002511124486 / Material number: 03810180 / Old Material Number: 68741301D

Product classification	
Thermal disinfectors	•
Machine affiliation	
PG 8581	•
PG 8591	•
Product category	
Insert/module for dental instruments	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of dental instruments	•
Capacity	
Hoses [number]	6
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	205
External dimensions, net width in mm	390
External dimensions, net depth in mm	450
External dimensions, gross height in mm	250
External dimensions, gross width in mm	485
External dimensions, gross depth in mm	495
Net weight in kg	1.40
Gross weight in kg	1.41