


A 13

Lid

For fixing delicate wash items.



- Lid in stainless steel
- 1 mm wire, 8 mm mesh size, 4 mm frames
- As replacement for E 103, E 104, E 105 and E 139 inserts

EAN: 4002511117204 / Material number: 03810200 / Old Material Number: 69745001D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 8505	•
Product category	
Lid	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	20
External dimensions, net width in mm	200
External dimensions, net depth in mm	320
External dimensions, gross height in mm	20
External dimensions, gross width in mm	220
External dimensions, gross depth in mm	300
Net weight in kg	0.14
Gross weight in kg	0.25