


UTS
Utensil box
for the ideal cleaning of small parts.



- Equipped with lid
- 3 compartments with a size of ?? x ??

EAN: 4002511117143 / Material number: 03810260 / Old Material Number: 69741001D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Other inserts	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	3
Area of application	
Reprocessing of MIS instruments	•
Reprocessing of surgical instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	93
External dimensions, net width in mm	102
External dimensions, net depth in mm	350
External dimensions, gross height in mm	110
External dimensions, gross width in mm	400
External dimensions, gross depth in mm	605
Net weight in kg	0.75
Gross weight in kg	1.00