

E 118

Insert

for the optimum loading of half petri dishes.



- Suitable for 38 half Petri dishes with Ø 100 mm
- equipped with 38 holders, height 70 mm, distance approx. 26 mm
- For use in upper basket

EAN: 4002511115767 / Material number: 03830270 / Old Material Number: 69511801D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 8505	•
<b>Product category</b>	
Insert/module for half petri dishes	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Half petri dishes [number]	38
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	120
External dimensions, net width in mm	460
External dimensions, net depth in mm	445
External dimensions, gross height in mm	160
External dimensions, gross width in mm	480
External dimensions, gross depth in mm	480
Net weight in kg	1.54
Gross weight in kg	1.90