


E 136

Insert

for optimum loading of 56 half petri dishes.



- Suitable for 56 half petri dishes with Ø 100 mm
- Equipped with 56 holders, at 70 mm height, spacing approx. 26 mm
- For use in lower basket
- The capacity can be doubled with insert E 137

EAN: 4002511115941 / Material number: 03830280 / Old Material Number: 69513601D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 8505	•
<b>Product category</b>	
Insert/module for half petri dishes	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
Distance between rails in mm	26
Rail height in mm	70
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Half petri dishes [number]	56
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	485
External dimensions, net width in mm	145
External dimensions, net depth in mm	445
External dimensions, gross height in mm	150
External dimensions, gross width in mm	450
External dimensions, gross depth in mm	535
Net weight in kg	2.10
Gross weight in kg	2.50