

E 143

Insert

For optimum accommodation of different instruments.



- 1 mm wire, 5 mm mesh size, 5 mm frames
- equipped with 2 hinged handles

EAN: 4002511115989 / Material number: 03830310 / Old Material Number: 69514301D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
PWD 8531	•
PWD 8532	•
PWD 8534	•
PWA 8672	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Mesh tray	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of dental instruments	•
Reprocessing of surgical instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	55
External dimensions, net width in mm	255
External dimensions, net depth in mm	230
External dimensions, gross height in mm	75
External dimensions, gross width in mm	285
External dimensions, gross depth in mm	260

E 143

Insert

For optimum accommodation of different instruments.

EAN: 4002511115989 / Material number: 03830310 / Old Material Number: 69514301D

Net weight in kg	0.49
------------------	------

Gross weight in kg	0.70
--------------------	------

