


A 2  
Cover net  
For fixing delicate wash items.



- Suitable for 1/2 inserts
- Plastic coated metal frame with plastic netting
- Size: 216 x 456 mm

EAN: 4002511115590 / Material number: 03830460 / Old Material Number: 69500201D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 8505	•
<b>Product category</b>	
Cover	•
<b>Product properties</b>	
Material	Wire, rilsan-coated, Plastic
Colour	Green, White, Black
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
Reprocessing of different utensils	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net width in mm	216
External dimensions, net depth in mm	456
External dimensions, gross height in mm	20
External dimensions, gross width in mm	230
External dimensions, gross depth in mm	465
Net weight in kg	0.38
Gross weight in kg	0.48