


A 3

Cover net

As cover for wash items in inserts.



- Plastic coated metal frame with plastic netting
- Suitable for 1/4 inserts
- 206 x 206 mm

EAN: 4002511115606 / Material number: 03830470 / Old Material Number: 69500301D

Product classification	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 6011	•
PLW 6111	•
PLW 8505	•
Product category	
Cover	•
Product properties	
Material	Wire, rilsan-coated, Plastic
Colour	Green, White, Black
Area of application	
Reprocessing of laboratory glassware	•
Reprocessing of different utensils	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	206
External dimensions, net depth in mm	206
External dimensions, gross height in mm	10
External dimensions, gross width in mm	210
External dimensions, gross depth in mm	200
Net weight in kg	0.30
Gross weight in kg	0.32