

AK 12
Insert
for holding different utensils.



- suitable for breast milk pumps, beakers & funnels, for example

EAN: 4002511115620 / Material number: 03830510 / Old Material Number: 69501201D

Product category	
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
PG 8562	•
Product category	
Insert/module for beakers	•
Insert/module for wide-necked glassware	•
Other inserts	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Reprocessing of different utensils	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	122
External dimensions, net width in mm	224
External dimensions, net depth in mm	434
External dimensions, gross height in mm	123
External dimensions, gross width in mm	315
External dimensions, gross depth in mm	505
Net weight in kg	0.74
Gross weight in kg	1.06