

E 442

Irrigation sleeves

for the optimum cleaning of MIS instruments.



- Suitable for instruments with Ø 4 – 8 mm, screw connection
- Length 121 mm mm, Ø 6 mm mm
- Cap, opening Ø 6 mm
- Mounting clip (mat. no. 4174850)

EAN: 4002511908208 / Material number: 04203830 /
Old Material Number: 69744201D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
PWA 8672	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Irrigation sleeve	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of MIS instruments	•
Reprocessing of lumened instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	121
External dimensions, gross height in mm	12
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	120
Diameter in mm	6
Net weight in kg	0.03
Gross weight in kg	0.03