

E 448

Silicone hose

for the optimum cleaning of various instruments.



- Luer lock adapter, male
- Length 300 mm mm, Ø 10 mm mm
- Jets, thread 8 x 1 mm

EAN: 4002511908475 / Material number: 04203880 / Old Material Number: 69744801D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Other accessories	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Transparent
Area of application	
Reprocessing of MIS instruments	•
Reprocessing of ophthalmology instruments	•
Reprocessing of lumened instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	300
External dimensions, net width in mm	8
External dimensions, net depth in mm	8
External dimensions, gross height in mm	10
External dimensions, gross width in mm	170
External dimensions, gross depth in mm	125
Diameter in mm	10
Net weight in kg	0.02
Gross weight in kg	0.02