

E 471
Spring clip
for injector jets with a Ø of 2.5 mm.



EAN: 4002512008464 / Material number: 04451930 / Old Material Number: 69747101D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Other accessories	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of anaesthetic instruments	•
Reprocessing of MIS instruments	•
Reprocessing of ENT instruments	•
Reprocessing of lumened instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, gross height in mm	12
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	120
Net weight in kg	0.00
Gross weight in kg	0.00