


E 803

Insert

for the optimum loading of up to 30 ear and nose specula.



- 160 compartments approx. 13 x 13 mm
- Mesh width on base: 1.7 mm
- For use in upper and lower basket

EAN: 4002512087735 / Material number: 04637750 / Old Material Number: 69580301D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8534	•
<b>Product category</b>	
Insert for ENT instruments	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	160
<b>Area of application</b>	
Reprocessing of ENT instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	63
External dimensions, net width in mm	165
External dimensions, net depth in mm	317
External dimensions, gross height in mm	80
External dimensions, gross width in mm	220
External dimensions, gross depth in mm	450
Net weight in kg	1.25
Gross weight in kg	1.30