

E 147/1

Insert

for the optimum loading of 10-12 mouthwash glasses.



- For 10-12 mouth rinse tumblers, max. Ø 80 mm
- Plastic-coated

EAN: 4002512122313 / Material number: 04746620 / Old Material Number: 68514703D

Product category	
Thermal disinfectors	•
Machine affiliation	
PG 8581	•
PG 8591	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Other inserts	•
Product properties	
Material	Wire, rilsan-coated
Colour	White
Area of application	
Reprocessing of dental instruments	•
Capacity	
Mouth rinsing cups [number]	12
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	155
External dimensions, net width in mm	220
External dimensions, net depth in mm	455
External dimensions, gross height in mm	260
External dimensions, gross width in mm	510
External dimensions, gross depth in mm	240
Net weight in kg	1.29
Gross weight in kg	1.59