

E 941
Mobile unit
for modules on 2 levels.



- Suitable for up to 2 modules per level depending on the size of the load
- Coupling of the mobile unit to the water and drying air supply
- The clearances for both levels are different

EAN: 4002512273411 / Material number: 04812530 / Old Material Number: 69594101D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Mobile unit for orthopaedic instruments	•
Other mobile units	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	421
External dimensions, net width in mm	619
External dimensions, net depth in mm	790
External dimensions, gross height in mm	700
External dimensions, gross width in mm	800
External dimensions, gross depth in mm	1200
Net weight in kg	16.30
Gross weight in kg	26.00
Companion pack	
Magnetic strip for automatic mobile unit recognition (incl. magnets) [number]	1