

E 467

Irrigation sleeves

For MIS instruments/clip forceps, length 205 mm, Ø 11 mm.



EAN: 4002512274340 / Material number: 05039640 /
Old Material Number: 69546701D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Irrigation sleeve	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of MIS instruments	•
Reprocessing of lumened instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	205
External dimensions, gross height in mm	60
External dimensions, gross width in mm	250
External dimensions, gross depth in mm	180
Diameter in mm	11
Net weight in kg	0.02
Gross weight in kg	0.02