

E 484
Insert
for holding different utensils



- Wire mesh 1.4 mm, mesh width 8 mm
- Can be equipped with holders: E 485, E 486, E 487, E 488, E 489

EAN: 4002512398718 / Material number: 05109190 / Old Material Number: 69548401D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Insert for kidney dishes/basins	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of theatre shoes	•
Reprocessing of basins	•
Reprocessing of kidney dishes	•
Reprocessing of surgical instruments	•
Reprocessing of different utensils	•
Capacity	
Theatre shoes [number]	16
Basins [number]	4
Kidney dishes [number]	9
Breathing masks [number]	9
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	150
External dimensions, net width in mm	470
External dimensions, net depth in mm	480
External dimensions, gross height in mm	560
External dimensions, gross width in mm	180
External dimensions, gross depth in mm	540
Net weight in kg	2.97
Gross weight in kg	3.78