


E 487

Bracket

for optimum loading of 4 theatre shoes up to a height of 280 mm.



- Several holders will fit in the E 484 insert

EAN: 4002512398749 / Material number: 05109250 / Old Material Number: 69548701D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Insert for theatre shoes	•
Area of application	
Reprocessing of theatre shoes	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Capacity	
Theatre shoes [number]	4
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	280
External dimensions, net width in mm	464
External dimensions, net depth in mm	10
External dimensions, gross height in mm	470
External dimensions, gross width in mm	470
External dimensions, gross depth in mm	25
Net weight in kg	0.37
Gross weight in kg	0.52