


E 488
Bracket
for the optimum loading of respiratory masks.



- Suitable for 3 respirator masks
- Several holders will fit in the E 484 insert

EAN: 4002512398756 / Material number: 05109260 / Old Material Number: 69548801D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Other inserts	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of different utensils	•
Capacity	
Breathing masks [number]	3
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	31
External dimensions, net width in mm	240
External dimensions, net depth in mm	210
External dimensions, gross height in mm	75
External dimensions, gross width in mm	280
External dimensions, gross depth in mm	510
Net weight in kg	0.27
Gross weight in kg	0.39