

## AUF 1 Holder

for optimum reprocessing of transmission instruments in upper basket O 177/1.



- Comprises: seal, insert nut, clamp
- Rack lower part, rack upper part, 5 filter plates

EAN: 4002512427517 / Material number: 05159510 / Old Material Number: 68551101D

<b>Product category</b>	
Thermal disinfectors	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8591	•
PWA 8672	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
<b>Product category</b>	
Holder and/or adapter for transmission instruments	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of dental instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, gross height in mm	125
External dimensions, gross width in mm	125
External dimensions, gross depth in mm	315
Net weight in kg	0.60
Gross weight in kg	0.94
<b>Companion pack</b>	
FP filter plate	•