

E 702

Attachment

to increase capacity by two further mesh trays.



- A further level for mobile unit E 525/2 or E 701/2
- For 2 further DIN mesh trays

EAN: 4002512400831 / Material number: 05221490 / Old Material Number: 69570201D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
Product category	
Other inserts/other modules	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Reprocessing of surgical instruments	•
Capacity	
DIN mesh trays [number]	2
SPRI mesh trays [number]	1
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	560
External dimensions, net width in mm	160
External dimensions, net depth in mm	530
External dimensions, gross height in mm	610
External dimensions, gross width in mm	260
External dimensions, gross depth in mm	535
Net weight in kg	3.50
Gross weight in kg	4.50