

E 496  
 Injector nozzle  
 with retaining clamp for intubation material.



EAN: 4002512365970 / Material number: 05235490 / Old Material Number: 69749601D

<b>Product category</b>	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
<b>Product category</b>	
Injector nozzle	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Area of application</b>	
Reprocessing of anaesthetic instruments	•
Reprocessing of MIS instruments	•
Reprocessing of lumened instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	120
External dimensions, gross height in mm	120
External dimensions, gross width in mm	30
External dimensions, gross depth in mm	30
Diameter in mm	4
Net weight in kg	0.75
Gross weight in kg	0.90