

E 496
 Injector nozzle
 with retaining clamp for intubation material.



EAN: 4002512365970 / Material number: 05235490 / Old Material Number: 69749601D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
Product category	
Injector nozzle	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of anaesthetic instruments	•
Reprocessing of MIS instruments	•
Reprocessing of lumened instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	120
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	120
External dimensions, gross width in mm	30
External dimensions, gross depth in mm	30
Diameter in mm	4
Net weight in kg	0.75
Gross weight in kg	0.90