

DOS 60/30

Dispenser pump

For dispensing liquid neutral/alkaline media.



- For dispensing liquid cleaning agent
- For dispensing chemical disinfectant
- For retro-fitting by Miele Service

EAN: 4002512417532 / Material number: 05267420 / Old Material Number: 69225003D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7825	•
G 7826	•
<b>Product category</b>	
Dispenser system	•
<b>Product properties</b>	
Production rate in ml/60 s	120
Length of suction lance in mm	333
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, gross height in mm	310
External dimensions, gross width in mm	250
External dimensions, gross depth in mm	430
Net weight in kg	1.70
Gross weight in kg	1.95