


ML  
Magnets  
for automatic recognition of mobile unit/spray arm monitoring



- 15 different configurable codes
- 6 magnets for the magnetic strip installed for each area of use

EAN: 4002512423083 / Material number: 05267440 / Old Material Number: 69745702D

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
Product category	
Magnetic strip for automatic mobile unit recognition	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel, Black
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	75
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	70
Net weight in kg	0.85
Gross weight in kg	1.00