

E 730
 Insert
 for the optimum loading of 20 theatre shoes.



EAN: 4002512423120 / Material number: 05271020 / Old Material Number: 69573001D

Product category	
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7825	•
G 7826	•
Product category	
Insert for theatre shoes	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of theatre shoes	•
Capacity	
Theatre shoes [number]	20
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	290
External dimensions, net width in mm	470
External dimensions, net depth in mm	540
External dimensions, gross height in mm	420
External dimensions, gross width in mm	540
External dimensions, gross depth in mm	600
Net weight in kg	8.70
Gross weight in kg	10.00