

AUF 2 Holder

for optimum reprocessing of transmission instruments in upper basket O 801/2.



- Comprises: seal, insert nut, clamp
- Rack lower part, rack upper part, 5 filter plates

EAN: 4002512427524 / Material number: 05282470 / Old Material Number: 68551701D

| | |
|---|-----------------|
| Product category | |
| Thermal disinfectors | • |
| Machine affiliation | |
| G 7831 | • |
| PG 8581 | • |
| PG 8591 | • |
| Product category | |
| Holder | • |
| Product properties | |
| Material | Stainless steel |
| Colour | Stainless Steel |
| Area of application | |
| Reprocessing of dental instruments | • |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0 |
| Frequency in Hz | 0 |
| Total rated load in kW | 0.00 |
| Fuse rating in A | 0 |
| Dimensions and weight | |
| External dimensions, net height in mm | 0 |
| External dimensions, net width in mm | 0 |
| External dimensions, net depth in mm | 0 |
| External dimensions, gross height in mm | 118 |
| External dimensions, gross width in mm | 122 |
| External dimensions, gross depth in mm | 305 |
| Net weight in kg | 0.21 |
| Gross weight in kg | 0.26 |
| Companion pack | |
| FP filter plate | • |