


O 801/2

Upper basket

for optimum cleaning with 10 nozzles and holding of various inserts.



- Front & back halves suitable for inserts
- Middle axis with 10 brackets and 10 nozzles Ø 4 mm, L 30 mm
- Vertical clearance: 200 mm
- With integrated spray arm

EAN: 4002512399463 / Material number: 05290520 / Old Material Number: 69180103D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
Product category	
Injector basket	•
Upper basket	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Maximum vertical clearance in mm	200
Area of application	
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Capacity	
Transmission instruments [number]	6
Lumened instruments [number]	10
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	267
External dimensions, net width in mm	381
External dimensions, net depth in mm	475
External dimensions, gross height in mm	480
External dimensions, gross width in mm	280
External dimensions, gross depth in mm	400
Net weight in kg	2.40
Gross weight in kg	3.10