

E 986  
Bracket  
for the optimum loading of 5 basins.



- For use in insert E 984

EAN: 4002512431347 / Material number: 05337770 / Old Material Number: 69598601D

<b>Product classification</b>	
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8527	•
PG 8528	•
<b>Product category</b>	
Insert for kidney dishes/basins	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of basins	•
<b>Capacity</b>	
Basins [number]	1
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	25
External dimensions, gross width in mm	330
External dimensions, gross depth in mm	435
Net weight in kg	0.19
Gross weight in kg	0.36