

E 520  
Insert  
for the optimum loading of root canal instruments.



- Suitable for 17 root canal instruments
- Hinged
- Can be steam sterilised at 121 °C/134 °C

EAN: 4002512472524 / Material number: 05409080 / Old Material Number: 68552001D

<b>Product category</b>	
Thermal disinfectors	•
<b>Machine affiliation</b>	
G 7831	•
PG 8581	•
PG 8591	•
<b>Product category</b>	
Other inserts	•
<b>Product properties</b>	
Material	Plastic
Colour	White
<b>Capacity</b>	
Root canal instruments [number]	17
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	45
External dimensions, net width in mm	75
External dimensions, net depth in mm	30
External dimensions, gross height in mm	80
External dimensions, gross width in mm	95
External dimensions, gross depth in mm	110
Net weight in kg	0.05
Gross weight in kg	0.10