


E 747

Module

for optimum loading of centrifuge tubes, phials and test tubes, for example.

EAN: 4002512591218 / Material number: 05464630 / Old Material Number: 69574701D



- 104 nozzles for centrifuge tubes, phials, test tubes, etc.
- 104 nozzles 2.5 x 110 mm

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
Product category	
Insert/module for test tubes	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of pipettes/phials	•
Capacity	
Test tubes [number]	104
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	168
External dimensions, net width in mm	492
External dimensions, net depth in mm	496
External dimensions, gross height in mm	650
External dimensions, gross width in mm	660
External dimensions, gross depth in mm	810
Net weight in kg	38.00
Gross weight in kg	44.00