


E 743

Module

for the optimum loading of narrow necked glassware 100 - 500 ml.



- With 36 nozzles (E 351) 4 x 160 mm with clips (E 353)

EAN: 4002512591164 / Material number: 05555250 / Old Material Number: 69574301D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
G 7825	•
G 7826	•
<b>Product category</b>	
Insert/module for narrow-necked glassware	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Narrow-necked glassware [number]	36
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	190
External dimensions, net width in mm	492
External dimensions, net depth in mm	496
External dimensions, gross height in mm	170
External dimensions, gross width in mm	530
External dimensions, gross depth in mm	560
Net weight in kg	9.70
Gross weight in kg	10.00