


E 746  
Module  
for the optimum loading of up to 23 pipettes.



- Suitable for 23 pipettes arranged diagonally
- 10 pipettes up to 560 mm and 13 pipettes up to 490 mm
- Distance between holders 20 or 26 mm

EAN: 4002512591195 / Material number: 05555280 / Old Material Number: 69574601D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
G 7825	•
<b>Product category</b>	
Insert/module for pipettes	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Capacity</b>	
Pipettes [number]	23
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	330
External dimensions, net width in mm	492
External dimensions, net depth in mm	496
External dimensions, gross height in mm	400
External dimensions, gross width in mm	590
External dimensions, gross depth in mm	540
Net weight in kg	24.60
Gross weight in kg	25.00