

DOS-P2  
Dispenser pump  
with hose and suction lance.



- Suitable for assembly in the dispenser pump casing DOS G2

EAN: 4002512615075 / Material number: 05573160 / Old Material Number: 59170080D

<b>Product category</b>	
Washing machines	•
<b>Machine/accessory category</b>	
PW 413	•
PW 413 Mop Star 130	•
PW 418	•
PW 5134 MOPSTAR 130	•
PW 6107	•
PW 6137	•
PW 5136	•
PW 6163	•
PW 6167	•
PW 6207	•
PW 6241	•
PW 6243	•
PW 6321	•
PW 6323	•
PW 811	•
PW 814	•
PW 818	•
PWM 500-08	•
PWM 500-09	•
PWM 506 Mop Star 60	•
PWM 507	•
PWM 508 Mop Star 80	•
PWM 509	•
PWM 900-09	•
PWM 906	•
PWM 907	•
PWM 908	•
PWM 511	•
PWM 909	•
PWM 514	•
PWM 520	•
PWM 912	•
PWM 916	•
PWM 920	•
PDW 909	•
<b>Product category</b>	
Dispenser system	•
<b>Product properties</b>	
Electrical connection	OHNE
<b>Dimensions and weight</b>	
External dimensions, net height in mm	540
External dimensions, net width in mm	540
External dimensions, net depth in mm	190
External dimensions, gross height in mm	210
External dimensions, gross width in mm	550
External dimensions, gross depth in mm	550
Net weight [kg]	1.30
Gross weight in kg	3.30

DOS-P2  
Dispenser pump  
with hose and suction lance.