


E 523

Insert

for the optimum loading of 6 mesh trays.



- Equipped with 7 holders (6 compartments), H 145, spacing 50 mm

EAN: 4002513199116 / Material number: 05649700 / Old Material Number: 69552301D

Product category	
Thermal disinfectors	•
Machine affiliation	
PG 8581	•
PG 8591	•
Product category	
Insert for trays	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	6
Area of application	
Reprocessing of dental instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	151
External dimensions, net width in mm	220
External dimensions, net depth in mm	450
External dimensions, gross height in mm	195
External dimensions, gross width in mm	240
External dimensions, gross depth in mm	495
Net weight in kg	1.40
Gross weight in kg	1.65