

E 470

Injector nozzle

for use in the mobile injector unit, Ø 2.5 x 90 mm.

EAN: 4002513248814 / Material number: 05701580 / Old Material Number: 69547001D



- Suitable for mobile injector units, with screw connection
- Suitable for bottles, Erlenmeyer flasks and round bottomed flasks for example

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Injector nozzle	•
Product properties	
Material	Wire, rilsan-coated, Stainless steel
Colour	Stainless Steel, White
Area of application	
Reprocessing of laboratory glassware	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	90
External dimensions, net width in mm	3
External dimensions, net depth in mm	3
External dimensions, gross height in mm	12
External dimensions, gross width in mm	100
External dimensions, gross depth in mm	150
Diameter in mm	3
Net weight in kg	0.02
Gross weight in kg	0.03