

E 470

Injector nozzle

for use in the mobile injector unit, Ø 2.5 x 90 mm.

EAN: 4002513248814 / Material number: 05701580 / Old Material Number: 69547001D



- Suitable for mobile injector units, with screw connection
- Suitable for bottles, Erlenmeyer flasks and round bottomed flasks for example

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
<b>Product category</b>	
Injector nozzle	•
<b>Product properties</b>	
Material	Wire, rilsan-coated, Stainless steel
Colour	Stainless Steel, White
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	90
External dimensions, net width in mm	3
External dimensions, net depth in mm	3
External dimensions, gross height in mm	12
External dimensions, gross width in mm	100
External dimensions, gross depth in mm	150
Diameter in mm	3
Net weight in kg	0.02
Gross weight in kg	0.03