


E 755

Module

for the optimum loading of narrow necked glassware 25 - 100 ml.



- Equipped with 36 nozzles (E 470) 2.5 x 90 mm with holder

EAN: 4002513248821 / Material number: 05701590 / Old Material Number: 69575501D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
G 7826	•
Product category	
Insert/module for narrow-necked glassware	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of pipettes/phials	•
Capacity	
Narrow-necked glassware [number]	36
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	130
External dimensions, net width in mm	492
External dimensions, net depth in mm	496
External dimensions, gross height in mm	545
External dimensions, gross width in mm	160
External dimensions, gross depth in mm	530
Net weight in kg	4.80
Gross weight in kg	5.40